





### Weva Face

#### Daily basic skin management system

- Can treat all people with different ages, skins, colors
- Can be used with any daily skin care products
- Can be used with any skin treat ment device






## Theory

### What is Weva Face?

Hydrafacial is a technology-driven rejuvenating skin treatment that delivers instant results you can see and feel. It promises to deliver long term skin health and can cater to the specific needs of all skin types.



It works by removing dead skin cells and enriching the skin with cleansing, hydrating and moisturising serums



It can be used for facial lift - ing, firming severe skin and diminishing wrinkles



# Handpiece

### Skin Peeling

Deep skin clean, removing skin cell detox, poretightening sebumcontrol & Whitehead.



### Lontophoresis



ION:Using (+) & (-) currents, gives skin elasticity with stimulating the skin tissue and muscles.

### Multipolar RF

Multipolar RF: Collagen remolding & improve skin elasticity .Safe and painless treatment with imme - diate effect.



### Ultrasound




Ultrasound: Melts fat, maintain skin texture smooth, absorb the nutritient deeply

# QUANTUM®

## Handpiece

### Cooling and Heating Handle

Cooling and heating system to contract pores and expand pores for helping solutions absorption



#### **BOHR Handle**

Use carbon oxygen bubble to release oxygen inside the skin for deep cleaning, remove dead skin cell to get smooth skin

### Oxygen Jet

Jet technology to infuse oxygen into the skin , anti-oxidation





# Applications

Deep cleaning

Skin Nourishment

Skin Refresh

Face Lifting


Skin rejuvenation

Facial skin care

Blackheads removal


Dead skin exfoliator removal







# Effect



Good Return		
Treatments	Procedures	Price
Facial deep cleaning	3-5	\$150
microdermabrasion	3-5	\$200
Blackhead Removal	3-5	\$200
Skin Brightening	3-5	\$150
Exfoliating	3-5	\$200
Skin Care Management	3-5	\$180
Shrinking Pores	3-5	\$160
Oxygen Renewal	3-5	\$200



# Technical Specification

Tornado tip	Large, Small
Vacuum range	Max 680 mmHg
Radio frequency	1Mhz, Multi polar(6Pol)
Lon lifting	500Hz (Digital Lon Lifting)
Ultrasound	1Mhz / 2W/cm <sup>2</sup>
Screen	10inch TFT
Input voltage	110V / 220V 50/60HZ
Noise level	45Db
Ultrasound working mode	Continous, Pulse
Net weight	39KG
Gross weight	71KG
Package Dimensions	59*53*135cm





### Testimonials



Dr. Erick Smith



My clients like the treatment very much. They enjoy the process very much. I am very happy to see my clients' skins get better!





Dra.Kelly Jones



"After the Hydra Facial' treatment, my patients notice that the skin is tighter, plumper and they are get ting compliments on their skin.





Dr.Albert Brown

Its cleaning ability is very good, and will not damage the skin. This machine should be a must-have for every beauty salon, it helps the fol low-up beauty project.

